



## MAHURKAR®\* Maxid™ Dual Lumen Chronic Care Catheter

### PATIENT BENEFITS

- Catheter design allows for navigation of difficult vasculature such as with femoral or left IJ placement.
- Designed for safety, performance and reliability
- Provides consistent flow rates that meet and exceed K/DOQI Guidelines of >300 mL/min.

### PERFORMANCE

- 14.5Fr Double-D\*\* lumen design optimizes flow with lower arterial pressures
- Carbothane™\* material provides additional kink resistance
- Silicone extensions provide kink resistance and allow view of blood
- Ultem™\* adapters resist cracking

### SAFETY

- Valved pull-apart safety sheath and VenaTrac™ over-the-wire insertion stylets reduce risk of air embolism and blood loss
- Safety tunneler allows for ease of insertion and reduces the risk of catheter detachment
- Large laser printed volumes are easy to read
- Clearly marked manufacturer's identification for future references





# MAHURKAR®\* Maxid™ Dual Lumen Chronic Care Catheter



MAHURKAR®\*  
Maxid Dual Lumen  
Catheter with  
sideholes



VenaTrac™ Over-the-  
Wire Insertion Stylets



Valved Pull-Apart  
Safety Sheath



Carbothane™\*  
material allows  
superior kink  
resistance, flexibility  
for placement, and  
maximum flows  
with low arterial and  
venous pressures.

|                                                                                                 |             | Priming Volumes      |               |          |   |
|-------------------------------------------------------------------------------------------------|-------------|----------------------|---------------|----------|---|
| Code                                                                                            | Description | Arterial (Red)       | Venous (Blue) | Qty Case |   |
| MAHURKAR®* Maxid™ Dual Lumen Catheter                                                           |             |                      |               |          |   |
| SPORT PACK                                                                                      | 8888119385  | 19/36 with sideholes | 1.6 mL        | 1.7 mL   | 5 |
|                                                                                                 | 8888123435  | 23/40 with sideholes | 1.8 mL        | 1.9 mL   | 5 |
|                                                                                                 | 8888128485  | 28/45 with sideholes | 2.0 mL        | 2.1 mL   | 5 |
|                                                                                                 | 8888133535  | 33/50 with sideholes | 2.2 mL        | 2.3 mL   | 5 |
| KITS                                                                                            | 8888145251  | 19/36 with sideholes | 1.6 mL        | 1.7 mL   | 5 |
|                                                                                                 | 8888145252  | 23/40 with sideholes | 1.8 mL        | 1.9 mL   | 5 |
|                                                                                                 | 8888145253  | 28/45 with sideholes | 2.0 mL        | 2.1 mL   | 5 |
|                                                                                                 | 8888145254  | 33/50 with sideholes | 2.2 mL        | 2.3 mL   | 5 |
| TRAYS                                                                                           | 8888145255  | 19/36 with sideholes | 1.6 mL        | 1.7 mL   | 1 |
|                                                                                                 | 8888145256  | 23/40 with sideholes | 1.8 mL        | 1.9 mL   | 1 |
|                                                                                                 | 8888145257  | 28/45 with sideholes | 2.0 mL        | 2.1 mL   | 1 |
|                                                                                                 | 8888145258  | 33/50 with sideholes | 2.2 mL        | 2.3 mL   | 1 |
| MAHURKAR®* Maxid™ Dual Lumen Catheter with VenaTrac™ Over-the-Wire Insertion Stylets (one pair) |             |                      |               |          |   |
| SPORT PACK                                                                                      | 8888119362  | 19/36 with sideholes | 1.6 mL        | 1.7 mL   | 5 |
|                                                                                                 | 8888123402  | 23/40 with sideholes | 1.8 mL        | 1.9 mL   | 5 |
|                                                                                                 | 8888128452  | 28/45 with sideholes | 2.0 mL        | 2.1 mL   | 5 |
|                                                                                                 | 8888133502  | 33/50 with sideholes | 2.2 mL        | 2.3 mL   | 5 |
| KITS                                                                                            | 8888145027  | 19/36 with sideholes | 1.6 mL        | 1.7 mL   | 5 |
|                                                                                                 | 8888145028  | 23/40 with sideholes | 1.8 mL        | 1.9 mL   | 5 |
|                                                                                                 | 8888145029  | 28/45 with sideholes | 2.0 mL        | 2.1 mL   | 5 |
|                                                                                                 | 8888145030  | 33/50 with sideholes | 2.2 mL        | 2.3 mL   | 5 |
| TRAYS                                                                                           | 8888145031  | 19/36 with sideholes | 1.6 mL        | 1.7 mL   | 1 |
|                                                                                                 | 8888145032  | 23/40 with sideholes | 1.8 mL        | 1.9 mL   | 1 |
|                                                                                                 | 8888145033  | 28/45 with sideholes | 2.0 mL        | 2.1 mL   | 1 |
|                                                                                                 | 8888145034  | 33/50 with sideholes | 2.2 mL        | 2.3 mL   | 1 |

| Components                           | Sport |      |          |
|--------------------------------------|-------|------|----------|
|                                      | Qty   | Pack | Kit Tray |
| 14.5Fr MAHURKAR®* Maxid™ Catheter    | 1     | •    | •        |
| 16Fr Valved Pull-Apart Safety Sheath | 1     | •    | •        |
| Tunneling Stylet, 21 cm              | 1     | •    | •        |
| 12Fr Tissue Dilator                  | 1     | •    | •        |
| 14Fr Tissue Dilator                  | 1     | •    | •        |
| Sealing Caps                         | 2     | •    | •        |
| 18G x 7 cm Introducer Needle         | 1     |      | •        |
| J/Straight .038" x 120 cm Guidewire  | 1     |      | •        |
| #11 Scalpel                          | 1     |      | •        |
| Telfa™ Island Dressings              | 2     |      | •        |
| 4 x 4 Cotton Gauze Sponges           | 4     |      | •        |
| 12 mL Syringe                        | 1     |      | •        |
| 6 mL Luer-Lock Syringes              | 2     |      | •        |
| 3 mL Luer-Lock Syringes              | 2     |      | •        |
| Needles (18G, 22G, 25G)              | 3     |      | •        |
| 2.0 Suture with Curved Needle        | 1     |      | •        |
| Fenestrated Drape                    | 1     |      | •        |
| 5 mL Ampule Lidocaine, 1%            | 2     |      | •        |
| Needle Holder                        | 1     |      | •        |
| Packet PVD Swabs                     | 1     |      | •        |
| PVD Ointment                         | 1     |      | •        |

®\* MAHURKAR and Double-D are registered trademarks of Sakharan D. Mahurkar, MD, used under license.

™\* Carbothane is a trademark of Noveon IP Holdings Corporation.

™\* Ultem is a trademark of Sabic Innovative Plastics IP B.V.